

CASTING TAPE

- Techform Casting Tape is engineered to produce an extremely strong, durable finished cast every time with great value, strength and durability.

Techform Standard - Product Highlights

- 100% fiberglass fabric ensures a strong and durable cast.
- Non-tack resin system and soft, conformable rolls allow for easy handling.
- Excellent lamination ensures a strong and durable cast.
- Heat treated fabric eliminates fraying when cut.
- Fantastic molding and conformability.
- Finger holes in the pouch allow for easy opening.



TECHFORM STANDARD

Specifications

Techform Standard

Fabric	100% fiberglass, heat-treated
Resin	Standard, no foam
Precision-welded ends	No
Fray-resistant cutting	Yes
Fabric Weight	Medium
Finished cast surface	Rough

Indication: - Used for any cast application



PART INFORMATION

WIDTH	WH	BK	OR	RD	PK	HP	LB	DB	PR	GN	RB
3.8 cm	154WH										
5 cm	204WH	204BK	204OR	204RD	204PK	204H P	204LB	204DB	204PR	204GN	204RB
7.6 cm	304WH	304BK	304OR	304RD	304PK	304H P	304LB	304DB	304PR	304GN	304RB
10 cm	404WH	404BK	404OR	404RD	404PK	404H P	404LB	404DB	404PR	404GN	404RB
12.7 cm	504WH										