

## CONTINUOUS PASSIVE MOTION (CPM)

### FLEX-MATE K500 KNEE CPM

- The Flex-Mate K500 is an anatomical Knee CPM that provides a range of motion (ROM) of minus 10 in extension to 120 of flexion.
- No configuration change is needed for right or left applications.
- At only 26 pounds (11.8kg), the Flex-Mate K500 is Lightweight and portable, making it ideal for both home and hospital use.
- The easy - to - use controls located on the handheld pendant allow the user to start or stop the unit, as well as adjust for flexion, extension, speed, and pause settings.
- In addition, a lockout feature can be incorporated to prevent unintentional or unauthorized setting changes.
- The Flex-Mate K500 accommodates a wide range of patient sizes due to its ability to be easily converted for both adult and pediatric applications with no additional accessories required.



PART#	DESCRIPTION
-------	-------------

06500	Flex-Mate K500 Knee CPM
-------	-------------------------

SOFT GOODS (PADDING)	
PART#	DESCRIPTION

06490	Flex-Mate K500 Pad Set
-------	------------------------

