

NIGHT SPLINT & AFO



With Slip-Resistant Tread

FORMFIT® NIGHT SPLINT

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon. Adjusts from 10°- 90° dorsiflexion.
- Patented 3-dimensional molding creates superior comfort and fit with padding where you need it.
- Soft, lycra-lined foam interior is plush, washable and breathable to eliminate heat build-up.
- Lycra-lined, foam-padded straps maintain comfort across the flexor tendons.
- User-friendly, center-release buckles simplify application.
- Lightweight, low-profile shell is sturdy and breathable.

PART#	SIZE	MEN	WOMEN	U/M
50023	Small	up to 7	up to 7.5	1 ea
50025	Medium	7.5 - 10	8 - 10.5	1 ea
50027	Large	10.5 & up	11 & up	1 ea

Indication: - Plantar Fasciitis
 - Achilles Tendonitis
 - Other lower extremity overuse injuries



AIRFORM® NIGHT SPLINT

- The Airform night splint provides an effective dorsal treatment for Plantar Fasciitis, Achilles Tendonitis, Drop Foot and Post-static pain.
- Its lightweight, low profile design makes it more comfortable. Better comfort means better compliance.
- Slips on and adjusts easily.
- Fits left or right foot, comes in 2 sizes.

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE	U/M
50045	Medium	37-43cm	37-43cm	1 ea
50047	Large	43.5-48.5cm	43.5-48.5cm	1 ea

Indication: - Plantar Fasciitis
 - Achilles Tendonitis
 - Drop Foot
 - Post-Static Pain



AFO LEAF SPRING

AFO Leaf Spring is a pre-fabricated polypropylene ankle-foot orthosis designed to support flaccid drop foot. It provides a semi-rigid section for toe clearance and support, yet has a thin, flexible foot plate for easy trimming. The absence of a heel section makes the Leaf Spring more comfortable to wear and provides a better fit in shoes.

- Injection-molded polypropylene is lightweight.
- Variable thickness throughout the orthosis provides strength.
- Provides good toe clearance and support.
- Thin, flexible foot part is easy to trim with a pair of scissor.
- Excellent fit for most types of shoes.
- Absence of heel section for improved comfort.

PART#	SIZE	HEIGH	LEG	FOOT LENGTH
O-110127	Small/Medium	27cm	Right	24cm
O-110227	Small/Medium	27cm	Left	24cm
O-110130	Large	30cm	Right	28cm
O-110230	Large	30cm	Left	28cm

Indication: - Drop foot secondary to CVA
 - Mild drop foot secondary to other neurological pathologies

NIGHT SPLINT & AFO

NAVIGAIT

NAVIGAIT™ is a new alternative for users with a mild foot drop. The amount of dorsiflexion is easily adjusted with an elastic coil to accommodate the user's specific needs.

- Developed to better assist user's with a mild foot drop to gain better dorsiflexion.
- The padded shell is positioned comfortably on the shin.
- Removable interface.
- The amount of dorsiflexion is adjusted through the position in the lacing and with the velcro attachment on the shin plate.
- Low profile design makes it invisible under slacks.
- NAVIGAIT™ is available in two sizes.
- The plastic shin part can be trimmed with scissors if needed.

Intended use and indications

NAVIGAIT™ is intended to support a foot when the ability to actively dorsiflex is reduced due to conditions such as Stroke, Multiple sclerosis, Post polio syndrome.

Contraindications

NAVIGAIT™ should not be used by patients who present with spasticity or severe oedema.

Limitations

NAVIGAIT™ does not address malpositions of the foot or foot deformities. These should be cared for with shoes and/or orthotics.

PART#	DESCRIPTION	SIZE	HEIGHT
28410 0011	NAVIGAIT™	Small/Medium	27.5 cm
28410 0013	NAVIGAIT™	Large/XLarge	31.5 cm



FOOT-UP®

Foot-up is a lightweight ankle-foot orthosis that offers dynamic support for drop foot or similar complaints. This revolutionary orthosis provides visible improvement in gait by supporting the foot the moment it is raised.

- Cushioned ankle wrap is made from a breathable, 3-layer material that can be worn comfortably for long periods of time.
- Plastic inlay fits discreetly between the tongue and laces of the shoe and is practically invisible when worn.
- Shoeless wrap allows to wear without shoes.
- Available in both Black and Beige colors.

ANKLE WRAPS

BLACK#	LEG	SIZE	ANKLE CIRCUMFERENCE
07810-1	Universal	Medium	18-21 cm
07810-2	Universal	Large	18-21 cm
07810-3	Universal	XLarge	27-33 cm

SHOELESS WRAPS

BLACK#	LEG	SIZE	ACCESSORY	FOOT CIRCUMFERENCE
27810-1	Universal	Small	Sholess Wrap	< 20 cm
27810-2	Universal	Medium	Sholess Wrap	20-24 cm
27810-3	Universal	Large	Sholess Wrap	24-28 cm
27810-4	Universal	XLarge	Sholess Wrap	> 28 cm

Sholess Wrap: measure the circumference of the foot immediately behind the ball of the foot.

Indication: - Dynamic support for drop foot or similar complaints for which Support of dorsiflexion is desirable.



* Measure the circumference of the ankle around the malleoli; measure the circumference of the foot immediately after the ball of the foot.