

SHOULDER

SHOULDER SUPPORTS - SlingShot 2

The SlingShot 2 features a comfortable, breathable Airmesh sling and a 15° abduction pillow. Quick release shoulder and waist strap buckles make this product easy to apply. Includes exercise ball to stimulate circulation and a thumb rest to minimize migration.

Features

- Comprised of Airmesh, a moisture wicking fabric, for enhanced breathability and comfort.
- Sizes (S, M & L)
- Universal left or right
- Quick release shoulder and waist strap

PART#	SIZE	FOREARM LENGTH
08502	Small	11.5" – 13" (29 – 33 cm)
08503	Medium	13.5" – 14.5" (34 – 37 cm)
08504	Large	15" – 16" (38 – 41 cm)

*Measurement taken from olecranon to knuckles

Indication: Rotator cuff repairs - Anterior repairs - Posterior dislocations
- Capsular shifts - Global shoulder instability



SHOULDER SUPPORTS - SlingShot 3

The SlingShot 3 is a clinician-driven innovation in shoulder bracing. Its design provide superior comfort for patient recovering from shoulder surgery while providing multiple options of post-operative support.

Features

- Innovative offloading shoulder harness for superior patient comfort.
- Convertible abduction pillow for post-operative options (15° abduction, 90° neutral and 45° neutral).
- Comprised entirely of Airmesh® fabric for enhanced breathability
- Lined with a moisture wicking blue fabric for maximum comfort
- Quick release shoulder and waist straps
- Sizes (S, M & L)
- Universal left or right
- Includes exercise ball

PART#	SIZE	FOREARM LENGTH
00042	Small	11.5" – 13" (29 – 33 cm)
00043	Medium	13.5" – 14.5" (34 – 37 cm)
00044	Large	15" – 16" (38 – 41 cm)

*Measurement taken from olecranon to knuckles

Indication: Rotator cuff repairs - Anterior repairs - Posterior dislocations
- Capsular shifts - Global shoulder instability



SHOULDER SUPPORTS - NEUTRAL WEDGE

Neutral Wedge is an immobilizer that positions the humerus in the neutral position. It includes a quick release immobilizing waist strap and comfortable neck pad.

Product Features

- Universal left or right
- Easy to fit pillow, sling and straps
- Size (S / M & L / XL)

PART#	SIZE	FOREARM LENGTH
11917	S / M	< 14.5" (37 cm)
11918	L / XL	> 14.5" (37 cm)

*Measurement taken from olecranon to knuckles

Indication: Bankart repairs - SLAP lesion repairs - Anterior dislocations
- Posterior capsule repairs - Anterior repairs - Rotator cuff repairs
- Total shoulder reconstruction



SHOULDER



SHOULDER SUPPORTS Shoulder Abduction Pillow

Shoulder Abduction Pillow is a shoulder immobilizer designed for varying degrees of abduction (10°-75°). When pillow is inflated, arm straps can be applied to limit posterior shift of the shoulder following rotator cuff repairs.

Features

- inflatable bladder for abduction of 10°-75°
- Limits posterior shift of the shoulder
- Universally sized abduction pillow
- Universal left or right

PART#	DESCRIPTION
01851	Shoulder Abduction Pillow (Universal)

Indication: Large rotator cuff repairs



SHOULDER ABDUCTION SLING

- Provides shoulder abduction at 10° with support
- Pillow portion rests on waist to support weight of arm
- Snap-on sling design for easy application
- Easy strap length adjustability
- Breathable cotton provides a comfortable fit
- Fits left and right

PART#	SIZE	FITS LENGTH*	U/M
204053	Small	Up to 31.7cm	1 ea
204055	Medium	31.7 - 39.3cm	1 ea
204057	Large	39.3cm-up	1 ea

* Measure elbow crease to fingertip

SHOULDER

BERGSCHULTZ UNIVERSAL ARM SLING

- Designed by Jim Bergschultz, R.O.T.
- Blue heavy denim material with contact closure
- One size only

PART#	SIZE
204070	Universal



UNIVERSAL ARM SLING - ADULT & CHILD

Provides general arm support and restricts movement to assist in healing sore or injured arms and elbows. Support left or right arm.

PART#	
204000-A	Universal Arm Sling -Adult
204000-C	Universal Arm Sling -Child



Indication: Shoulder injuries or Instabilities

