

WALKER

REBOUND® AIR WALKER

- Full circumferential shell with additional semi-soft toe bumper for frontal toe protection.
- Soft Flex Edge® reduces pressure points.
- Full rocker bottom is especially engineered to help promote a natural gait.
- Stretch-to-fit shell for optimal fit.
- Shell can be trimmed for additional stretch in the calf.
- Moisture-wicking softgood helps keep skin dry and comfortable.
- Linked medial/lateral bladders for easy inflation and equal compression.
- Integrated pump and air release valve for easy access.



- 1 Fully-integrated air pump and release valve provides easy, tool-free inflation and deflation for fine-tuned, consistent compression.
- 2 Flex Edge™ overmolding conforms to each patient's unique anatomy to help limit the formation of pressure points.
- 3 Rocker bottom promotes a natural, stable gait, while a trimmable, semi-soft toe bumper provides added protection.
- 4 Ventilated panels wick moisture to prevent microbial buildup and improve patient comfort.



Standard Version



Low Top Version

STANDARD VERSION

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE
B-242900001	XSmall	up to 35.5cm	up to 35.5cm
B-242900002	Small	36-39cm	36-38.5cm
B-242900003	Medium	40-44.5cm	39-43cm
B-242900004	Large	44.5-46.5cm	44.5-46.5cm
B-242900005	XLarge	46.5-48.5cm	46cm +

LOW TOP VERSION

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE
B-242900061	XSmall	up to 35.5cm	up to 35.5cm
B-242900062	Small	36-39cm	36-38.5cm
B-242900063	Medium	40-44.5cm	39-43cm
B-242900064	Large	44.5-46.5cm	44.5-46.5cm
B-242900065	XLarge	46.5-48.5cm	46cm +

HEEL WEDGE*

PART#	SIZE	DESCRIPTION
B-242902501	XSmall - XLarge	Heel Wedge

- Indication:
- Grade 2 and 3 ankle sprains
 - Stable fractures
 - Post-operative
 - Post-trauma
 - Rehabilitation
 - Soft tissue injuries



HEEL WEDGE*

- Used to offload the Achilles Tendon
- Maintains the angle of plantar flexion at 5°, 10°, 15°, 20° or 25°

WALKER

EQUALIZER® PREMIUM AIR WALKER

- Patented pneumatic system increases stability while decreasing pain and swelling.
- Air is easily adjusted to achieve desired compression and custom fit.
- Accommodates the swelling patterns that occur throughout the rehabilitation process.

STANDARD VERSION

PART BLACK#	SIZE	MEN'S SHOE	WOMEN'S SHOE
A-W0200BLK	XSmall	2-4	3.5-5.5
A-W0400BLK	Small	4.5-7	6-8
A-W0600BLK	Medium	7.5-10.5	8.5-11.5
A-W0800BLK	Large	10.5-12.5	11.5-13.5
A-W1000BLK	XLarge	12.5+	13.5+

- Indication:**
- Soft tissue injuries, grade 2 and 3 ankle sprains
 - Stable fractures
 - Post-operative
 - Trauma
 - Rehab use



Standard Version

EQUALIZER® LOW TOP WALKER

LOW TOP VERSION - BLACK

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE
W0300BLK	XSmall	2-4	3.5-5.5
W0500BLK	Small	4.5-7	6-8
W0700BLK	Medium	7.5-10.5	8.5-11.5
W0900BLK	Large	10.5-12.5	11.5-13.5
W1100BLK	XLarge	12.5+	13.5+

- Indication:**
- Soft tissue injuries, grade 2 and 3 ankle sprains
 - Stable fractures
 - Post-operative
 - Trauma
 - Rehab use



Low Top Version

RANGE OF MOTION WALKER

- The Range of Motion Joint can be fixed in any position between 0-45° of plantar flexion.
- The range of motion joint adjusts in 7.5° increments in both dorsi and plantar flexion, and can either be completely locked out at any desired interval, or left with desired range.
- Available in Black only.

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE
RW0400	Small	4.5 – 7	6 – 8
RW0600	Medium	7.5 – 10.5	8.5 – 11.5
RW0800	Large	10.5– 12.5	11.5 – 13.5

- Indication:** - Ideal for Achilles tendon repair

